

DRILLED- END CARDIOVASCULAR TAPER CUTTING NEEDLES

EXAMPLE OF NEEDLE CODE

VHSYZ 180 20 51 28 DRS

1 TYPE

V - Drilled end,
Cardiovascular

2 SHAPE

D - 3/8 Circle
H - 1/2 Circle
X - Compound curve

3 POINT

SY - Round body
taper cutting

5 ANGLE OF CURVATURE

135° - (3/8 Circle)
160° - (1/2 Circle)
165° - (1/2 Circle)
170° - (1/2 Circle)
175° - (1/2 Circle)
180° - (1/2 Circle)

6 LENGTH

Length of needle
(in mm)

8	24
10	25
12	26
13	28
15	30
18	35
20	36
22	

7 WIRE SIZE

Wire dia. of needle in
hundredths of mm.

20	51
23	56
28	58
30	61
33	66
38	71
40	81
43	86
46	100
48	

9 GRADE OF STEEL

D - AISI 302

11 SILICONISED

B - Black
S - Siliconised

8 HOLE SIZE

Diameter of drilled hole
in hundredths of mm.

9	30
10	32
11	33
12	36
13	40
14	42
15	43
16	44
17	46
19	47
20	48
21	49
23	51
25	53
27	55
28	

4 SPECIAL FEATURE (IF REQD.)

(VIEW D)

Z - Square body

10 RIBBED BODY

R - Ribbed body

