

DRILLED- END STERNAL TAPER CUTTING NEEDLES

EXAMPLE OF NEEDLE CODE

WHSY165 49 157 95 ARS

1 TYPE

W - Sternal needles

2 SHAPE

H - 1/2 Circle

3 POINT

SY - Taper cutting

4 ANGLE OF CURVATURE

165° - (1/2 Circle)
180° - (1/2 Circle)

5 LENGTH

Length of needle
(in mm)

48
49
50
58
59
60

6 WIRE SIZE

Wire dia. of needle in
hundredths of mm.

130
140
157

8 GRADE OF STEEL

A - AISI 420

7 HOLE SIZE

Diameter of drilled hole
in hundredths of mm.

63	93
65	95
73	103
84	105
85	

9 RIBBED BODY

R - Ribbed body

10 SILICONISED

S - Siliconised

